

Fast, Accurate, Easy to Use

Gluco Spark

Blood Glucose Measuring System



Features:


- 0.5 μ L Blood Sample Volume.
- 5 Seconds Fast Test Time.
- No Coding.
- Large LCD Display.
- Strip Ejection.
- Date and Time

CE IVD

ISO 15197 - 2013 ISO 9001 - 2015 ISO 13485 - 2016

SENSA  CORE

Strip Advantages



Specifications

Glucose ranges	:	20 - 800 mg/dL
Analysis time	:	5 sec
Sample volume	:	0.5 microliters
Sample	:	Whole blood
Sample mode	:	Capillary
Interference	:	Minimal
Fill detection	:	Automatic
Display	:	48 mm (Diagonally) LCD
Shelf life of strips	:	24 months (from manuf. date)
Strips packing	:	25, 50 and 100 (2 Vials x 50 Strips)
Memory	:	500 Readings
Strips Storage Temp.	:	0-50°C
Operating Temp.	:	2-50°C
Relative Humidity	:	10-95% (Non Condensing)
Battery	:	3v Lithium Battery (CR2032), Lasts for 1000 tests

Soft Inserting Test Strip

Inserts the test strip easily and smoothly.

Gluco Spark is india's first, next generation, highly advanced ultra low volume glucose monitoring system available at affordable price. This glucometer uses highly advanced biosensor technology and requires 0.5 microliters sample volume and presents result in 5 sec.

Gluco Spark uses carbon electrode with laser patterning to ensure precise results. **Gluco Spark** strips are manufactured with tight tolerances and do not require any coding. **Sensa Core's** biosensor technology solves a number of immediate technical problems related to glucose strip by significantly enabling sensitive and accurate performance of glucose strip in a cost-effective way.

Manufactured by :

Sensa Core Medical Instrumentation Pvt. Ltd.
Plot No: 3, E P I P Zone, Pashamylaram, Sangareddy (Dist.), Hyderabad,
Telangana - 502307, INDIA. Tel: +91-8455-223400, www.sensacore.com

SENSA CORE
Healthcare and Diagnostic Products