

Certificate of Analysis

Certificate No. SC/QA/SOP/004/01

Release Date: 01.02.2017

Product Name : Daily Cleaner

Product Code : ST-EA-CN-001

Storage Conditions : 2-8 ° C

Quantity : 15 ml x 6 alkaline aqueous solution

Function : To Remove protein deposited on flow tube of Electrolyte Analyzer.


Quality Specifications : Passed.

Handling : Exercise normal laboratory precautions.

Stability : Package stability (expiration date) is listed on the product label.

Reagent Contents : Alkaline Solution contains preservative in water.
For *in vitro* diagnostic use only.

Certified By,
Dr. Kalachar HCB (Ph.D)


Quality Assurance
Chemistry division

Sensa Core Medical Instrumentation Pvt. Ltd.

Plot, No: 3, Export Promotion Industrial Park, Pashamylaram, Sangareddy (Dist), 502307, Telengana, INDIA.

Tel: +08455 - 223400, www.sensacore.com