

Certificate of Analysis

Certificate No. SC/QA/SOP/004/02

Release Date:28.04.2016

Product Name : Quality controls (QC Level-1, 2 and 3 and Urine Level 1 & 2)

Product Code : ST-EA-CN-002

Storage Conditions : 18-25 ° C

Quantity : 12 × 10 ml each (QC 1 - 3 No's., QC 2 - 4 No's, QC 3 - 3 No's, Urine Level 1 and Level 2 one bottle each)

Function : For monitoring the performance of Sensa Core Electrolyte Analyzer.

Quality Specifications : Passed.

Handling : Exercise normal laboratory precautions.

Stability : Package stability (expiration date) is listed on the product label.

Reagent Contents : Buffered solution of electrolytes (Sodium, Potassium, Chloride, Lithium, Calcium). This control contains no human or biological materials. For *in vitro* diagnostic use only.

Certified By,



Dr. Kalachar HCB (Ph.D)
Quality Assurance
Chemistry division

Sensa Core Medical Instrumentation Pvt. Ltd.

Plot, No: 3, Export Promotion Industrial Park, Pashamylaram, Sangareddy (Dist), 502307, Telengana INDIA.

Tel: +08455 - 223400, www.sensacore.com