

Daily Cleaners Kit

For Use With

Electrolyte Analyzer Only.

Intended Use

For *in vitro* diagnostic use.

(Reagent-4) Daily Cleaner Diluent Alkaline Solution
6 bottles (6x15mL)

Storage and Stability

Store at 18-25⁰c until expiration date on labels

Irritant! Irritating to eyes and skin. Avoid contact with skin and eyes. In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.



For a definition of symbols used in product labeling, please refer to *Understanding the Symbols* in the Appendix section of the Operator's Manual.

Daily Cleaners Kit

For Use With

Electrolyte Analyzer Only.

Intended Use

For *in vitro* diagnostic use.

(Reagent-4) Daily Cleaner Diluent Alkaline Solution
6 bottles (6x15mL)

Storage and Stability

Store at 18-25⁰c until expiration date on labels

Irritant! Irritating to eyes and skin. Avoid contact with skin and eyes. In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.



For a definition of symbols used in product labeling, please refer to *Understanding the Symbols* in the Appendix section of the Operator's Manual.